

Adjustable Arm Rest



### Coburg Micro Surgery Stool

Suitable for all types of micro surgery including: Neuro Surgery, Ophthalmic Surgery & ENT Surgery.

Code: **3014/2016**

- Arm rests with swivel ball joints, detachable and autoclavable
- Hydraulic height adjustment with foot pump
- Adjustable back rest
- Foot brake to give stability
- Comfortable Surgeons seat, either round or saddle shaped
- 5 anti-static castors
- Height range 61 cm to 76cm



Available with round seat



### Coberg Stool

Suitable for all types of micro surgery including: Neuro Surgery, Ophthalmic Surgery & ENT Surgery.

Code: **3014**

- Hydraulic height adjustment with foot pump
- Adjustable back rest
- Foot brake to give stability
- Comfortable Surgeons seat, either round or saddle shaped
- 5 anti-static castors
- Height range 61 cm to 76cm



Available with round seat