


MICRO DISCECTOMY RETRACTOR SET



- This Specialist Retractor allows the Consultant the optimal exposure to the Disc with the minimum incision.
- The Design of the Retractor allows unhindered access around the wound with the jointed arms folding along the patients Back.
- The Black finish eliminates distracting reflections under the Microscope during Surgery.
- The wide selection of Blades available allows the Surgeon full flexibility in approaches to the Disc.

Basic Micro-Discectomy Retractor Set

SP5200-00

Complete set consisting of:

Qty	Order Code	Description
1 each	SP5202-00	Retractor Frame
4 pairs	SP5204-40	Wide Muscle Blades Sizes 4cm to 7cm depth
4 pairs	SP5206-40	Narrow Muscle Blades Sizes 4cm to 7cm depth
4 pairs	SP5208-40	Hook Blades Sizes 3cm to 6cm depth
4 pairs	SP5210-40	Toothed Blades Sizes 3cm to 6cm depth
1 each	SP5200-10	Instrument Container

Extended Micro-Discectomy Retractor Set

SP5200-01

Complete set consisting of:

Qty	Order Code	Description
2 each	SP5202-00	Retractor Frames
6 pairs	SP5204-60	Wide Muscle Blades Sizes 3cm to 8cm depth
6 pairs	SP5206-60	Narrow Muscle Blades Sizes 3cm to 8cm depth
6 pairs	SP5208-60	Hook Blades Sizes 2cm to 7cm depth
6 pairs	SP5210-60	Toothed Blades Sizes 2cm to 7cm depth
1 each	SP5200-10	Instrument Container

MICRO DISCECTOMY RETRACTOR SET



Retractor
Frame
with Jointed
arms

SP5202-00

Wide Muscle Blades

Dimensions



SP5204-03	3cm depth
SP5204-04	4cm depth
SP5204-05	5cm depth
SP5204-06	6cm depth
SP5204-07	7cm depth
SP5204-08	8cm depth

Hook Blades

Dimensions



SP5208-02	2cm depth
SP5208-03	3cm depth
SP5208-04	4cm depth
SP5208-05	5cm depth
SP5208-06	6cm depth
SP5208-07	7cm depth

Narrow Muscle Blades

Dimensions



SP5206-03	3cm depth
SP5206-04	4cm depth
SP5206-05	5cm depth
SP5206-06	6cm depth
SP5206-07	7cm depth
SP5206-08	8cm depth

Toothed Blades

Dimensions



SP5210-02	2cm depth
SP5210-03	3cm depth
SP5210-04	4cm depth
SP5210-05	5cm depth
SP5210-06	6cm depth
SP5210-07	7cm depth