

### “A” or “V” Pattern?

Throughout the series, the “V” pattern is the all-purpose clamp, commonly used on all veins, and almost all arteries. The “A” pattern, available in all sizes, is designed for exceptionally thick-walled vessels and slippery arteries. The slight incurve of the tips of the “A” pattern jaw provides extra security against slippage.

### The Best has become even Better!

What started as a laboratory clamp over thirty years ago, has become the most popular micro clamp due to the many advantages this clamp series has for its users:















- Light and compact design
- Extremely durable
- Corrosion resistant metallurgy
- A unique gripping surface
- Smooth sliding bar action on all approximator clamps
- Clamp profile with parallel sided jaws
- Choice of matt or black finish



Double Clamps

Single Clamps

### Ordering details

SINGLE CLAMPS		DOUBLE APPROXIMATING CLAMPS		DOUBLE CLAMPS WITH FRAME	
00396 V 00397 A	 CLAMP SIZE <b>B - 1</b> 0.4 - 1.0mm	00414 V 00415 A	 CLAMP SIZE <b>ABB - 1 1</b> 0.4 - 1.0mm	00408 V 00409 A	 CLAMP SIZE <b>ABB - 1</b> 0.4 - 1.0mm
00398 V 00399 A	 CLAMP SIZE <b>B - 2</b> 0.6 - 1.5mm	00416 V 00417 A	 CLAMP SIZE <b>ABB - 2 2</b> 0.6 - 1.5mm	00410 V 00411 A	 CLAMP SIZE <b>ABB - 2</b> 0.6 - 1.5mm
00400 V 00401 A	 CLAMP SIZE <b>B - 3</b> 1.0 - 2.25mm	00418 V 00419 A	 CLAMP SIZE <b>ABB - 3 3</b> 1.0 - 2.25mm	00412 V 00413 A	 CLAMP SIZE <b>ABB - 3</b> 1.0 - 2.25mm
00325 V	 CLAMP SIZE <b>HD - S</b> 1.5 - 3.5mm	00326 V	 CLAMP SIZE <b>HD - D</b> 1.5 - 3.5mm	 <p>Clamp Applicator CAF-4</p> <p>Clamp Box CB-1</p>	
00286	 CLAMP SIZE <b>RD - S</b> 2.0 - 5.0mm	00287	 CLAMP SIZE <b>RD - D</b> 2.0 - 5.0mm		