

CASPAR MICRO LUMBAR DISCECTOMY RETRACTOR SET



Caspar Micro Discectomy Retractor

- This Retractor Set allows a micro surgery approach to Lumbar Disc
- Specula have a black surface finish to prevent glare
- Fenestrated blades on the specula reduce the risk of movement during surgery

Complete Set

SP5100-00

Consisting of:

Lateral Blades	Depth	Specula	Depth
SP5101-40	40mm	SP5100-40	40mm
SP5101-45	45mm	SP5100-45	45mm
SP5101-50	50mm	SP5100-50	50mm
SP5101-55	55mm	SP5100-55	55mm
SP5101-60	60mm	SP5100-60	60mm
SP5101-65	65mm	SP5100-65	65mm
SP5101-75	75mm	SP5100-75	75mm
SP5101-80	80mm	SP5100-80	80mm
SP5101-85	85mm	SP5100-85	85mm
SP5101-10	100mm	SP5100-10	100mm

} Optional extra deep

- SP5102-00 Counter Retractor Frame
- SP5103-00 Blade Ejector Forceps
- SP5104-00 Sterilising Container
- SP5105-00 Tray Insert for Specula
- SP5106-00 Tray Insert for Lateral Blades

CASPAR MICRO LUMBAR DISCECTOMY RETRACTOR SET

Lateral Blades



Blade Ejector Forceps

SP5103-00



Specula Range



Caspar Micro Discectomy Retractor Frame

SP5102-00

