

“A” or “V” Pattern?

Throughout the series, the “V” pattern is the all-purpose clamp, commonly used on all veins, and almost all arteries. The “A” pattern, available in all sizes, is designed for exceptionally thick-walled vessels and slippery arteries. The slight incurve of the tips of the “A” pattern jaw provides extra security against slippage.















- Light and compact design
- Extremely durable
- A unique gripping surface
- Smooth sliding bar action on all approximator clamps
- Clamp profile with parallel sided jaws
- Matt or black finish



Acland
Clamp Container
MV4400-00



Ordering details

SINGLE CLAMPS	DOUBLE APPROXIMATING CLAMPS	DOUBLE CLAMPS WITH FRAME
00396 V 00397 A  CLAMP SIZE B-1 0.4 - 1.0mm	00414 V 00415 A  CLAMP SIZE ABB - 11 0.4 - 1.0mm	00408 V 00409 A  CLAMP SIZE ABB-1 0.4 - 1.0mm
00398 V 00399 A  B-2 0.6 - 1.5mm	00416 V 00417 A  ABB - 22 0.6 - 1.5mm	00410 V 00411 A  ABB - 2 0.6 - 1.5mm
00400 V 00401 A  B-3 1.0 - 2.25mm	00418 V 00419 A  ABB - 33 1.0 - 2.25mm	00412 V 00413 A  ABB - 3 1.0 - 2.25mm
00325  HD-S 1.5 - 3.5mm	00326  HD - D 1.5 - 3.5mm	 <p>Clamp Box CB-1</p> <p>Clamp Applicator CAF-4</p>
00286  RD-S 2.0 - 5.0mm	00287  RD - D 2.0 - 5.0mm	

