



## Caspar Vertebral Body Distraction Set

Complete set consisting of:

### SP3001-00

- For retraction of the vertebral body in the cervical spine to allow access to the disc
- Less instrumentation in the disc space
- Distraction screws are secured into the vertebral body either side of the disc, the retractor is then slid over the screws and retracted

Qty	Order Code	Description
1	SP3002-00	Retractor Right
1	SP3003-00	Drill Guide Right
1	SP3004-00	Screw Driver
1	SP3005-00	Distraction Drill
1	SP3006-14	Sterile Distraction Screws 14mm Boxes of 10

## Distraction Retractor

SP3002-00 Right (illustrated)

SP3002-10 Left



## Drill Guide

SP3003-00 Right (illustrated)

SP3003-10 Left



## Screw Driver

SP3004-00



## Distraction Drill

SP3005-00



- For pre-drilling pilot holes for screws with drill guide 2mm

## Distraction Screws

- 7cm post length
- Sterile packed



	Thread length
SP3006-12	12mm
SP3006-14	14mm
SP3006-16	16mm



Supplied Sterile in boxes of 10