

OPERATING CHAIRS



- Variable height adjustment
- Comfortable seat and back supports
- Arm rests & foot brake for extra stability during micro surgery

MERCIAN
Surgical Instruments of Excellence

Adjustable Arm Rest



Coburg Micro Surgery Stool

Suitable for all types of micro surgery including: Neuro Surgery, Ophthalmic Surgery & ENT Surgery.

3014/2016

- Arm rests with swivel ball joints, detachable and autoclavable
- Hydraulic height adjustment with foot pump
- Adjustable back rest
- Foot brake to give stability
- Comfortable Surgeons seat, either round or saddle shaped
- 5 anti-static castors
- Height range 61 cm to 76 cm



Available with round seat



Coberg Stool

Suitable for all types of micro surgery including: Neuro Surgery, Ophthalmic Surgery & ENT Surgery.

3014

- Hydraulic height adjustment with foot pump
- Adjustable back rest
- Foot brake to give stability
- Comfortable Surgeons seat, either round or saddle shaped
- 5 anti-static castors
- Height range 61 cm to 76 cm



Available with round seat

Back rest



**The Severn Surgeons
Operating Stool with Back Rest**
Suitable for all types of surgery

222.167

- Comfortable Surgeons seat
- Adjustable back rest
- 5 anti-static castors
- Height adjustment with foot plate height range 55cm to 78cm

Foot plate for
height adjustment

**The Severn Surgeons
Urology Operating Stool
with Back Rest**

222.157

- Lower height range 48cm to 67cm



**The Avon Surgeons Operating
Stool without Back Rest**
Suitable for all types of surgery

222.165

- Comfortable Surgeons seat
- 5 anti-static castors
- Hydraulic height adjustment with foot plate, range 55cm to 78cm

Foot plate for
height adjustment



**The Avon Surgeons
Urology Operating Stool
with Back Rest**

222.155

- Lower height range 48cm to 67cm





MERCIAN

Surgical Instruments of Excellence

established over
40 years

Mercian Surgical Supply Co. Limited
10 Topaz Business Park, Topaz Way, Bromsgrove, Worcestershire B61 0GD. UK.
Tel: +44 (0) 844 879 1133 Fax: +44 (0) 844 879 1155 Email: info@merciansurgical.com
www.merciansurgical.com