


BREAST RETRACTORS

Tebbetts Fibre Optic Breast Retractor

With fibre optic illumination and suction



- Fibre optic illumination with suction for smoke evacuation
- Angled Retractor tip with blunt teeth
- Polished mirror finish on blade to reflect light inside the pocket

With Teeth

Code	Length	Width
260-109-240	90mm	24mm
260-109-300	90mm	30mm
260-136-150	150mm	36mm
260-136-190	190mm	36mm

Blunt

Code	Length	Width
260-190-240	90mm	24mm
260-190-300	90mm	30mm
260-243-150	150mm	36mm
260-136-019	190mm	36mm

FIBRE OPTIC LIGHT LEAD

Fibre Optic Light Lead

Code	Length	Diameter
250-035-180	1800mm	3.5mm Ø
250-035-230	2300mm	3.5mm Ø



Adaptors for Light Lead

Retractor Fitting



Lightsource Fitting

