

LUMBAR SPINAL SET


LEVEL 3



LEVEL 2



LEVEL 1



Mercian introduce a new set of spinal instruments for lumbar spine surgery.

The Mercian Lumbar Spine Set has three levels of instruments consisting of Bone Currettes and Nerve Root Retractors on the first level, Bone Nibblers and Hooks on the second level with the lower level consisting of Kerrison Punches 2mm to 5mm and Disc Ronguers.



SP1000-00 Complete Set

Complete set consisting of:

Qty	Order Code	Description	
Titanium Kerrison Punches - x 4			
1	FSKF 4082 ETI	Kerrison Titanium Punch, 2mm - 20cm/8" length	LEVEL 1
1	FSKF 4083 ETI	Kerrison Titanium Punch, 3mm - 20cm/8" length	
1	FSKF 4084 ETI	Kerrison Titanium Punch, 4mm - 20cm/8" length	
1	FSKF 4085 ETI	Kerrison Titanium Punch, 5mm - 20cm/8" length	
Disc Ronguers x 6			
1	C7-S	Cushing Disc Ronguer, 2mm - 18cm/7" straight	LEVEL 1
1	LG7-S	Love-Gruenwald Disc Ronguer, 3mm - 18cm/7" straight	
1	S7-S	Spurling Disc Ronguer, 4mm - 18cm/7" straight	
1	LG7-U	Love-Gruenwald Disc Ronguer, 3mm - 18cm/7" upward	
1	S7-U	Spurling Disc Ronguer, 4mm - 18cm/7" upward	
1	SP1110-23	Williams Disc Ronguer, toothed 1mm - 23cm/9" straight	
Bone Nibblers x 6			
1	SP1018-24	Leksell Bone Nibbler, 3mm bite size - 24cm/9½" length	LEVEL 2
1	SP1019-24	Leksell Bone Nibbler, 5mm bite size - 24cm/9½" length	
1	SP1003-18	Jansen Bone Ronguer D/A, curved 5mm bite size - 18cm/7" length	
1	SP1027-21	Zielke Disc Ronguer, 5mm bite size - 21cm/8¼" length	
1	SP1030-21	Rat Toothed Nibbler, 8mm bite size - 21cm/8¼" length	
1	SP1006-22	Stille-Luer Bone Nibbler, 10mm bite size - 22cm/8½" length	
Bipolar Forceps			
1	70-01-92	Bayonet Bipolar Forceps, 1.5mm tips x 23cm/9" length	LEVEL 2
Ball Tipped Probes x 2			
1	SP5250-13	Ball End Probe with 4mm tip - 23cm/9" length	LEVEL 2
1	SP5250-14	Ball End Probe with 6mm tip - 23cm/9" length	
Bone Curettes x 8			
1	SP8040-27	Curette, straight, cup size 0 (2mm) - 27cm/10" length	LEVEL 3
1	SP8050-27	Curette, straight, cup size 1 (3mm) - 27cm/10" length	
1	SP8140-27	Curette, angled up, cup size 0 (2mm) - 27cm/10" length	
1	SP8150-27	Curette, angled up, cup size 1 (3mm) - 27cm/10" length	
1	SP8250-23	Bruns Curette, straight, cup size 2 (4mm) - 23cm/9" length	
1	SP8240-23	Bruns Curette, straight, cup size 1 (5mm) - 23cm/9" length	
1	SP8350-23	Bruns Curette, angled up, cup size 2 (4mm) - 23cm/9" length	
1	SP8340-23	Bruns Curette, angled up, cup size 1 (5mm) - 23cm/9" length	
Cobb Elevators x 2			
1	SP7561-13	Cobbs Elevator, 13mm width x 32cm/12¾" length	LEVEL 3
1	SP7561-19	Cobbs Elevator, 19mm width x 32cm/12¾" length	
Adson Hooks x 2			
1	SP9101-20	Adson Sharp Hook 90°, 5mm tip - 20cm/8" length	LEVEL 3
1	SP9102-20	Adson Blunt Hook 90°, 5mm tip - 20cm/8" length	
Dissectors x 2			
1	SP9206-23	MacDonald Dissector - 23cm/9" length	LEVEL 3
1	SP9203-23	Watson Cheyne Dissector - 23cm/9" length	
Nerve Root Retractors x 4			
1	SP5041-21	Krayenbuhl Nerve Root Retractor - 21cm/8¼" length	LEVEL 3
1	SP5042-24	Krayenbuhl Nerve Root Retractor - 24cm/9½" length	
1	SP5250-02	McCulloch Nerve Root Retractor 90° (Level 2)	
1	SP5013-31	Nerve Root Retractor with Suction - 31cm/12" length	