



Heavy Orthopaedic Mallet

- 25cm/10" length
- Head size 40mm Ø
- Weight 900g/32oz
- Ergonomic handle

HA5806-21



A very heavy mallet suitable for hip surgery with an ergonomic handle for comfort

Orthopaedic Mallet

- 23cm/9" length
- **Regular** head size 35mm Ø
- Weight **550g/20oz**

HA5804-23

- **Large** head size 40mm Ø
- Weight **700g/25oz**

HA5805-23



Offered in two weights suitable for large bone surgery

Heath Mallet

- 19cm/7½" length
- Stainless steel head
- Head size 40mm Ø
- Weight 1020g/36oz

HA5809-19



Traditional mallet for hip surgery

Ombredanne Mallet

- 24cm/9½" length
- Head size 40mm Ø
- Weight 550g/20oz

HA5806-24



Stainless Steel Mallet

- 22cm/8½" length
- Head size 35mm Ø
- Weight 500g/18oz
- 1 nylon face

HA5813-22



Nylon face is replaceable

Toffee Mallet

- 19cm/7½" length
- Head size 22mm Ø
- Weight 165g/6oz
- Lead filled head
- Solid stainless steel head

HA5808-21



Ideal for small bone surgery hand or facial surgery

Mercian Surgical Supply Co. Limited

10 Topaz Business Park, Topaz Way, Bromsgrove, Worcestershire B61 0GD. UK.

www.merciansurgical.com

 @MercianSurgical  Mercian Surgical

Tel: +44 (0) 844 879 1133 Fax: +44 (0) 844 879 1155 Email: info@merciansurgical.com