

# ORTHOPAEDIC K-WIRES, TWIST DRILLS & STEINMANN PINS



**MERCIAN**

Surgical Instruments of Excellence


[www.merciansurgical.com](http://www.merciansurgical.com)

# TWIST DRILLS

## Round Fitting Twist Drill

- Stainless steel
- Supplied sterile in boxes of 10


Code	Diameter	Length	Flute
TD1045-11	1.1mm Ø	45mm	20mm
TD1070-15	1.5mm Ø	70mm	30mm
TD1085-20	2.0mm Ø	85mm	38mm
TD1185-20	2.0mm Ø	185mm	38mm
TD1095-25	2.5mm Ø	95mm	43mm
TD1195-25	2.5mm Ø	195mm	50mm
TD1130-32	3.2mm Ø	130mm	40mm
TD1095-35	3.5mm Ø	95mm	42mm
TD1195-35	3.5mm Ø	195mm	50mm
TD1130-45	4.5mm Ø	130mm	55mm



## AO Fitting Twist Drill

- Stainless steel
- Supplied sterile in boxes of 10

Code	Diameter	Length	Flute
TD2060-11	1.1mm Ø	60mm	12mm
TD2085-15	1.5mm Ø	85mm	18mm
TD2100-20	2.0mm Ø	100mm	20mm
TD2110-20	2.0mm Ø	110mm	32mm
TD2110-25	2.5mm Ø	110mm	32mm
TD2110-25G	2.5mm Ø	110mm	32mm (Gold tip)
TD2145-25	2.5mm Ø	145mm	43mm
TD2100-27	2.7mm Ø	100mm	30mm
TD2150-27	2.7mm Ø	150mm	50mm
TD2145-32	3.2mm Ø	145mm	43mm
TD2195-32	3.2mm Ø	195mm	50mm
TD2110-35	3.5mm Ø	110mm	40mm
TD2195-35	3.5mm Ø	195mm	50mm
TD2145-45	4.5mm Ø	145mm	50mm
TD2195-45	4.5mm Ø	195mm	50mm



- Stainless steel
- Supplied sterile in boxes of 10




Single Ended	100mm/4"	125mm/5"	150mm/6"
Ø 0.6mm	KW6100 - 06	KW6125 - 06	KW6150 - 06
Ø 0.7mm	KW6100 - 07	KW6125 - 07	KW6150 - 07
Ø 0.8mm	KW6100 - 08	KW6125 - 08	KW6150 - 08
Ø 0.9mm	KW6100 - 09	KW6125 - 09	KW6150 - 09
Ø 1.0mm	KW6100 - 10	KW6125 - 10	KW6150 - 10
Ø 1.1mm	KW6100 - 11	KW6125 - 11	KW6150 - 11
Ø 1.2mm	KW6100 - 12	KW6125 - 12	KW6150 - 12
Ø 1.4mm	KW6100 - 14	KW6125 - 14	KW6150 - 14
Ø 1.5mm	KW6100 - 15	KW6125 - 15	KW6150 - 15
Ø 1.6mm	KW6100 - 16	KW6125 - 16	KW6150 - 16
Ø 1.8mm	KW6100 - 18	KW6125 - 18	KW6150 - 18
Ø 2.0mm	KW6100 - 20	KW6125 - 20	KW6150 - 20



Double Ended	100mm/4"	125mm/5"	150mm/6"
Ø 0.6mm	KW7100 - 06	KW7125 - 06	KW7150 - 06
Ø 0.7mm	KW7100 - 07	KW7125 - 07	KW7150 - 07
Ø 0.8mm	KW7100 - 08	KW7125 - 08	KW7150 - 08
Ø 0.9mm	KW7100 - 09	KW7125 - 09	KW7150 - 09
Ø 1.0mm	KW7100 - 10	KW7125 - 10	KW7150 - 10
Ø 1.1mm	KW7100 - 11	KW7125 - 11	KW7150 - 11
Ø 1.2mm	KW7100 - 12	KW7125 - 12	KW7150 - 12
Ø 1.3mm	KW7100 - 13	KW7125 - 13	KW7150 - 13
Ø 1.4mm	KW7100 - 14	KW7125 - 14	KW7150 - 14
Ø 1.5mm	KW7100 - 15	KW7125 - 15	KW7150 - 15
Ø 1.6mm	KW7100 - 16	KW7125 - 16	KW7150 - 16
Ø 1.8mm	KW7100 - 18	KW7125 - 18	KW7150 - 18
Ø 2.0mm	KW7100 - 20	KW7125 - 20	KW7150 - 20

## Steinmann Threaded Pins

- Stainless steel
- Trocar point
- Fully threaded
- Supplied sterile in boxes of 10




	20cm/8"	23cm/9"
Single Ended		
Ø 1.4mm	KW8200-14FT	KW8230-14FT
Ø 1.6mm	KW8200-16FT	KW8230-16FT
Ø 1.8mm	KW8200-18FT	KW8230-18FT
Ø 2.0mm	KW8200-20FT	KW8230-20FT
Ø 2.4mm	KW8200-24FT	KW8230-24FT
Ø 2.8mm	KW8200-28FT	KW8230-28FT
Ø 3.2mm	KW8200-32FT	KW8230-32FT
Ø 4.8mm	KW8200-48FT	KW8230-48FT

	20cm/8"	23cm/9"
Double Ended		
Ø 1.4mm	KW9200-14FT	KW9230-14FT
Ø 1.6mm	KW9200-16FT	KW9230-16FT
Ø 1.8mm	KW9200-18FT	KW9230-18FT
Ø 2.0mm	KW9200-20FT	KW9230-20FT
Ø 2.4mm	KW9200-24FT	KW9230-24FT
Ø 2.8mm	KW9200-28FT	KW9230-28FT
Ø 3.2mm	KW9200-32FT	KW9230-32FT
Ø 4.8mm	KW9200-48FT	KW9230-48FT

## Denham Mid-Threaded Pins

- Stainless steel
- Trocar point/single ended
- Centrally threaded
- Supplied sterile in boxes of 5



	20cm/8"	23cm/9"
Single Ended		
Ø 4.0mm	KW8200-40CT	KW8230-40CT
Ø 4.5mm	KW8200-45CT	KW8230-45CT
Ø 5.0mm	KW8200-50CT	KW8230-50CT

Mercian Surgical Supply Co. Limited

10 Topaz Business Park, Topaz Way, Bromsgrove, Worcestershire B61 0GD. UK.

[www.merciansurgical.com](http://www.merciansurgical.com)

@MercianSurgical Mercian Surgical

Tel: +44 (0) 844 879 1133 Fax: +44 (0) 844 879 1155 Email: [info@merciansurgical.com](mailto:info@merciansurgical.com)